




Instructions for Use

Pro-Flex® LP Align

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	23
NO	Bruksanvisning	27
DA	Brugsanvisning	32
SV	Bruksanvisning	36
EL	Οδηγίες Χρήσης	41
FI	Käyttöohjeet	46
NL	Gebruiksaanwijzing	50
PT	Instruções de Utilização	55
PL	Instrukcja użytkowania	60
TR	Kullanım Talimatları	64
RU	Инструкция по использованию	69
JA	取扱説明書	74
ZH	中文说明书	78
KO	사용 설명서	82
CS	Návod k použití	86
SK	Indikácie použitia	91
RO	Indicații de utilizare	95
HR	Upute za uporabu	100
HU	Használati javallat	104
BG	Указания за употреба	109
SL	Indikácie pre použitie	114
UA	Інструкція з використання	118

DESCRIPTION

The device is a prosthetic foot with user adjustable heel height. It is composed of an integrated hydraulic ankle unit, with a male pyramid and three composite blades. It features a sandal toe and a full-length composite blade. The bottom blade is secured with heel bolts.

The device consists of the following components (**Fig. 1**):

- A. Male Pyramid
- B. Lock/Release Button
- C. Upper Housings
- D. Top Blade
- E. Middle Blade
- F. Bottom Blade Bolts
- G. Bottom Blade
- H. Lower Housing
- I. Posterior Spacer
- J. Anterior Spacer Wedge

INTENDED USE

The device is intended as a part of a prosthetic system that replaces the foot and ankle function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to moderate impact use, e.g., walking.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: The device must NEVER be used when Lock/Release Button is unlocked.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

DEVICE SELECTION

Verify that selected variant of the device is suitable for the impact level and weight limit according to the following table.

Warning: Do not exceed weight limit. Risk of device failure.

Incorrect category selection may also result in poor device function.

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	221-256
Low Impact Level							
Size 22-24	1	1	2	3	4	5	N/A

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	221-256
Low Impact Level							
Size 25-27	1	1	2	3	4	5	6
Size 28	3	3	3	3	4	5	6
Moderate Impact Level							
Size 22-24	1	2	3	4	5	N/A	N/A
Size 25-27	1	2	3	4	5	6	N/A
Size 28	3	3	3	4	5	6	N/A

ASSEMBLY INSTRUCTIONS

Foot Cover and Spectra Sock (Fig. 2)

Caution: To avoid pinching fingers, always use a shoehorn.

1. Put the foot into the Spectra Sock.
2. Remove the foot cover attachment.
3. Use the applicable (straight) end of a shoehorn to put the foot with the Spectra Sock into the Foot Cover.
4. Move the shoehorn up to fully push the foot into the Foot Cover.
5. Put the foot cover attachment back onto the Foot Cover.
6. After alignment is complete, fix the Spectra Sock to the prosthesis to seal against dust and dirt.

Note: The Spectra Sock must be pulled up to prevent it from interfering with moving parts of the foot.

If required doff the Foot Cover as follows:

1. Remove the foot cover attachment.
2. Insert the applicable (edged) end of a shoehorn behind the foot.
3. Push the shoehorn down and pull the foot out of the Foot Cover.
4. Fully remove the Spectra Sock.

Prosthesis

Assemble prosthesis with applicable devices.

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

Warning: Ensure proper attachment by following the applicable device assembly instructions. Remove the protective film on the pyramid after fitting.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 3)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at either the patellar tendon level or at the ischial tuberosity level (D).
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect either the socket or the knee to the foot and establish the correct knee center height.
3. If using a prosthetic knee: Position knee according to knee alignment instructions

4. On the lateral side of the socket, make a first mark (alignment reference point) at the midpoint of the socket at either patellar tendon level or the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through the alignment reference point (D).
6. Set the correct socket angles for flexion/extension and abduction/adduction.
7. If using a prosthetic knee: Use the applicable adapters to connect the socket to the prosthesis.

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Ensure the patient is familiar with the functioning of the device.

The heel to toe action can be influenced by:

- Heel Stiffness
- Anterior-Posterior positioning of device
- Dorsi-Plantarflexion
- Shoe Characteristics

Consider the following actions if needed:

Symptoms

- Device comes to flat position too early (patient feels he/she is sinking into a hole)
- Climbing over the toe requires extra energy
- Toe feels too stiff
- Knee hyperextends

Action

- Check Wedge instructions (if applicable)
- Shift socket anterior (or device posterior)
- Consider dorsiflexion
- Check heel of the shoe and shoe performance

Symptoms

- Rapid heel to toe movement
- Poor control over prosthesis at initial contact
- Minimal energy return feeling
- Too little push off from the toe
- Knee becomes unstable

Action

- Check Wedge instructions (if applicable)
- Shift socket posterior (or device anterior)
- Consider plantarflexion
- Check heel of the shoe and shoe performance

USAGE

Lock/Release Mechanism

- Press the medial (inner) side of the Lock/Release Button to unlock the ankle.
Note: The device CANNOT be used in the unlocked position for weightbearing.
- Adjust the ankle to the appropriate heel height by pressing on the toe or heel (it is also adjustable by hand when the prosthesis is not fitted, and resting on a table).
- The pylon or shin section should be perpendicular to the floor when the foot is in the shoe.
- When appropriate heel height adjustment is achieved, press the lateral (outer) side of the Lock/Release Button to lock the ankle.

- Ensure the Lock/Release Button is locked prior to weight bearing.

Note: The Lock/Release Button may be stiff at first, or when not used for long periods of time. Pushing in and out several times will loosen the button. Similarly, the ankle may also be stiff, and should be moved up and down several times to loosen.

Note: An extension for the Lock/Release Button should not be used, as it may result in accidental release.

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

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The Posterior and Anterior Spacers may be replaced if they become worn.

To replace:

- Remove the Foot Cover and Spectra Sock from the device.
- Compress the Posterior Spacer to expose the Anterior Spacer.*
- Remove the Anterior Spacer.
- Position the foot in maximum dorsiflexion to expose and remove the Posterior Spacer.
- Replace with new Posterior Spacer, ensuring the spacer is fully seated.
- Compress the Posterior Spacer to expose the Anterior Spacer.*
- Replace with new Anterior Spacer, ensuring the spacer is fully seated.

*Place the wedge included in the replacement kit between the Bottom and Middle Blades to ease Posterior Spacer compression.

Note: Ensure the wedge is removed before using the device.

Noise from Foot

Noise may occur if sand or debris is present in device. In that case, the healthcare professional should doff the foot, clean it with the help of compressed air and replace the Spectra sock if it is damaged.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to two million load cycles. Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 Label		
Category	Weight (Kg) Moderate Impact	Label Text
1	52	ISO 10328-P3*-52 kg
2	59	ISO 10328-P3*-59 kg
3	68	ISO 10328-P3-68 kg
4	77	ISO 10328-P3-77 kg
5	88	ISO 10328-P4-88 kg
6	100	ISO 10328-P5-100 kg
*Actual test loads reflect maximum body mass		

ISO 10328 - "P" - "m"kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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