

Descrierea produsului

Balance Foot S is a lightweight, waterproof foot specifically designed to meet the needs of low active walkers. For these users, the cushioning heel foam with gradual stiffening and wide sole blade provides enhanced stability perception, and the perception of a smooth rollover. The anatomical foot cover features a steady-grip sole and a sandal clamping toe, for the opportunity to wear different types of footwear.

Össur recommends using Balance Foot S in conjunction with the sleeveless elevated vacuum suspension system Unity®, which adds minimal weight and no build height, and does not interfere with Balance Foot S dynamics. Unity is designed to provide secure suspension in combination with Iceross Seal-In® technology, while stabilizing volume fluctuations and maintaining comfortable socket fit throughout the day.

General information

- Level of limb difference:

Dezarticulație de genunchi, Dezarticularea șoldului, Transfemural and Transtibial

- Impact Level:

Scăzut

- Characteristics

Pârghie completă până la vârf, Multi-axial ground compliance, Vârf tip sanda, Antepicior divizat and Rezistent la apă

- Available sizes

22-30

- Patient weight range

45 – 147 kg (99 – 324 lbs)

- Max patient weight

147 kg (324 lbs)

Characteristics

Pârghie completă până la vârf

Pârghia completă până la vârful lamei/degetului se potrivește cu lungimea piciorului sănătos, asigurând un mers mai lin și simetric.

Multi-axial ground compliance

Multi-axial compliance of the foot provides improved ground contact during ambulation on uneven surfaces.

Vârf tip sanda

Înveliș îngust pentru picior, corect din punct de vedere anatomic, pentru utilizare cu o încălțăminte cu vârf tip sandală.

Antepicior divizat

Caracteristica de antepicior divizat permite labei piciorului să se adapteze la suprafața de dedesubt, persoanele cu amputații putându-se deplasa natural pe o suprafață neuniformă.

Rezistent la apă

Fie că este vorba de apă dulce, apă sărată sau apă clorurată, produsele Össur rezistente la apă vă vor oferi protecția maximă de care aveți nevoie. Utilizarea materialelor de primă clasă oferă protecție permanentă împotriva pătrunderii nocive a apei, chiar și după scufundare, continuând să ofere funcționalitatea la care vă așteptați.